



## **Once notified each Employee has 24 hours to complete the test.**

Please remind your employees that they are required to bring picture identification with them to the testing facility. At the testing site they should identify themselves as part of the Screensafe testing pool. They will also need to retain the testing receipt the facility gives them which will need to be returned to you, the employer, to provide proof that the employee has complied with the testing request.

## **IMPORTANT! WE HAVE SOMETHING NEW TO OFFER, REGISTER ONLINE AND SAVE SOME TIME**

**Please remember to be at the testing site at least 1 hour before closing.**

Many collection facilities now require a Barcode.

The Barcode will eliminate confusion regarding who the person is testing for and why.

1. Register online: The participant can register themselves online  
It will take approximately two (2) minutes to schedule.

Go to [www.escreengo.com](http://www.escreengo.com)

Enter Scheduling Code:

**752SSI**

The participant will receive a text message or email with the Barcode.

The participant need only show that Barcode and a valid photo ID. to the collection facility staff.

2. Contact ScreenSafe to set up a location call tollfree 877-727-3369, we will go through the same process, and the member will receive a barcode via text message or email.

3. Some locations still allow for walk-in service, go to

<https://www.screensafeinc.com/page/search-sites>

and type in the zip code. Location information and instructions will be listed.

LAB CORP collection site please say testing for ScreenSafe and give the account number:

**LabCorp Account 235582**